




SENSUM



HELIX

Automatic Lifting
and Dedusting
of Tablets & Capsules

Computer
Vision
Systems

HELIX

Automatic Lifting and Dedusting of Tablets & Capsules

Sensum | 02





Sensum HELIX is a general-purpose lifting and dedusting machine. HELIX is compact, mobile and modular so that it can be easily transferred and adapted to various production stages:

- tablet pressing
- capsule filling
- tablet coating
- capsule weighting
- quality inspection
- blister filling
- bottle filling
- and other ...








Wherever and whenever there is need for gentle and efficient lifting and/or dust removing from tablets or capsules.

Highlights

- Double function: lifting & dedusting
- Product independent: tablets & capsules
- Gentle product manipulation
- Compact and mobile
- Input container of arbitrary size
- Modular construction
- Simple cleaning and lifting to arbitrary height
- Stand alone or integrated with SPINE inspection machine



Technical data

	Application	• Tablets and capsules
	Speed	• 1.000.000 products/hour
	Electrical connection	• 230V 1/N/PE 2A
	Compressed air	• 4-6 Bar, 60 Nm ³ /h
	Dimensions (LxWxH)	• ~ 630 x 700 x 760 - 2740 mm
	Weight	• ~ 170 kg
	Construction	• CE, GMP, FDA acceptable materials



SENSUM

Sensum, Computer
Vision Systems
Tehnološki park 21
SI - 1000 Ljubljana
Slovenia, EU
www.sensum.eu
info@sensum.eu



Specialised in high quality assurance in pharmaceutical industry by integrated computer vision systems.

R&D collaboration with the Imaging Technologies Lab, Department of Electrical Engineering, University of Ljubljana.

30 years of experience in image processing and analysis.